




24/7 PUBLIC ACCESS DEFIBRILLATORS in MACDUFF

Banff and Macduff
Heartstart, PAD and
Community First Responders

E: bandm.heartstart.cfr.foch@gmail.com
and on Facebook

for use in case of Cardiac Arrest



1. Call 999. Start CPR. Send someone to the nearest PAD;
2. Pull break glass on green box to get code; alarm will sound;
3. Enter code on keypad and open cabinet door;
4. Take Defibrillator to the patient. Turn on and follow instructions.